

Please cite the J.B. Lawrence Lab or publications for use of this protocol.  
Thanks!

---

### **Antibody Staining + RNA Hybridization**

Usually we do the antibody reaction first, but some antibodies only work if they come last (same protocol but reversed).

#### **Primary AB**

1. Rinse coverslips in 1xPBS 10 min
2. Make appropriate dilution of antibody in 1xPBS/1% BSA and add 1unit/ul of RNAsin. Place a 50-80µl drop onto a glass plate lined with parafilm
3. Place coverslip, cell side down, on top of drop
4. Cover with a 2<sup>nd</sup> piece of parafilm, seal the sides like an envelope to prevent evaporation and incubate 1 hour at 37°C.
5. Wash coverslips at room temp:
  - 1xPBS – 10min on shaker
  - 1xPBS + 0.1% Triton – 10min on shaker
  - 1xPBS – 10min on shaker

#### **Secondary AB**

6. Make a appropriate dilution of conjugated ( Fitc, Texas Red, Alexa 488 or 594) antibody in 1XPBS/ 1% BSA add 1unit/ul of RNAsin.
7. Place drops of secondary antibody onto a parafilm lined glass plate
8. Place coverslip on top of drop, cells down
9. Cover with a 2<sup>nd</sup> piece of parafilm, wrap entire plate with tin foil, (to keep reaction in the dark) and incubate 1 hour at 37°C.
10. Wash coverslips at room temp :
  - 1xPBS – 10min on shaker

1xPBS + 0.1% Triton – 10min on shaker  
1xPBS – 10min on shaker

11. Fix coverslips in 4% Paraformaldehyde for 10 min, room temp, in dark
12. Rinse twice in 1xPBS

### **RNA hybridization**

13. Follow Lawrence Lab standard RNA hybridization protocol from this website.
14. Use a detector with a different fluorochrome than the one used with the antibody.

### **DAPI DNA stain**

15. Incubate in DAPI stain, 30sec-1 min, in dark
16. Rinse twice with 1xPBS
17. Mount coverslips onto slides using Vectashield (Vector Labs) mounting media and seal edges with fingernail polish.
18. Slides are stored in a slide folder at -20°C